






How to Create a Recording

1. Turn on the Braille + by hitting the **X** or also known as the cancel button.
2. To do this, locate the space bar and move your finger up slightly to locate the **X** (or the third button under the down arrow ( **MENU X**)
3. Turn the Braille + vertically so the microphone is located at the top. The microphone is tactually identifiable as it is ridged.
4. To locate the record button move your hand up the left side of the Braille+ starting from the bottom. The first button contacted is the record button.
5. Press the record button and hold for about 2-3 seconds.
6. The device will say begin recording.
7. This is when you starting speaking your recorded message.
8. To pause while recording press the select key  located in the middle of the arrow keys.
9. To restart and continue recording press the select key  again.
10. The Braille + will say recording.
11. When you're finished recording press the ok button.  This is the first button located under the arrow keys (ok, menu, cancel) or ( **MENU X**)
12. The Braille + will state recording and the current name of the file, this is usually the current date and time, followed by saved.

To play back your new recording press the select key  located in the middle of the arrow keys.