







## How to Rename a Recording

1. Press the menu button located directly below the down arrow key.  
(ok, menu, cancel) or (  **MENU**  )
2. The Braille + will say options menu 1 rename recording.
3. Press the number 1 or the ok button.  This is the first button located under the arrow keys. (ok, menu, cancel) (  **MENU**  )
4. The Braille + will say "New name computer Braille required and current name of the file which is mostly like the date." This mean type the name using letters only not contractions. When using numbers a number sign is not required however the number must be in Nemeth.
5. Press the ok button  when you done typing the new name of the recording.
6. The Braille + will then announce the name of the file .mp3 and the position of the recording i.e. 1 of 7.